


RESTORATIONS WITH GLASS FIBER POSTS


GF Posts & Precision Drills

Fiberglass posts with high strength and reliability



SEM image of sectioned Polydentia GF Posts. Glass fibers are well distributed and with a high density (80% of total surface section).

Image by Laboratory of Biophysics, University of Lille - France



Polydentia GF Posts
Translucent, optimal light transmission

Clinical pictures by courtesy of Dr. G. Derchi - La Spezia, Italy



Image 1
Initial

Image 2
Clinical application

Image 3
Final

ADVANTAGES

- Biocompatible
- An alternative to metal posts
- Excellent and strong adhesion to composites
- For highest aesthetic results
- Easy to remove in case of revision

DESCRIPTION

Polydentia GF Posts consist of a large number of unidirectional, individual fiberglasses (diameter: 20 µm), embedded in a polyester resin matrix - a material especially appreciated for its allergy safety. The parallel, longitudinal fibers reinforce the post's structure and enhance its flexibility. The white, opaque color of Polydentia GF Posts ensures excellent aesthetic results, without undesirable coloring in the cervical area of the tooth. If endodontic retreatment becomes necessary, the post can be easily removed by drilling. The clinical relevant product properties are the following:

Endo-Design: parallel shape with tapered end for optimal adaptation to the anatomy of the root.

Complete range of sizes: the introkit contains a complete size offer, covering all possible posterior and anterior treatments. The full length is 20 mm (.78").

Composition: uniform distribution and high density of glassfibers (80% of total surface).


Radiopacity: tests based on ISO 4049 showed a result of 150% Al.

Translucency: high light transmission into the root canal.

Elasticity: no stress concentrations in dentine or root fractures due to a elasticity modulus of 42 GPa.

Mechanical properties: increased stability of the post-and-core restorations due to the high flexural strength (1500 MPa) and high mechanical resistance to compression (650 MPa).

Retention: active surface of the post, for optimal micro-mechanical adhesion, providing an excellent and strong adhesion to the core composite.



Description and content	REF	Description and content	REF
Polydentia GF Posts & Precision Drills Introkit 20 Assorted Posts, 5 each (ø 1.1 mm, 1.3 mm, 1.5 mm, 1.7 mm) + 4 Matching Precision Drills	5280	Polydentia GF Posts & Precision Drills Other Assortments 15 Assorted Posts, 5 each (ø 1.3 mm, 1.5 mm, 1.7 mm) + 3 Matching Precision Drills	5292
		15 Assorted Posts, 5 each (ø 1.1 mm, 1.3 mm, 1.5 mm) + 3 Matching Precision Drills	5285
		Polydentia GF Posts Other Assortments 15 Assorted Posts, 5 each (ø 1.3 mm, 1.5 mm, 1.7 mm)	5291
		15 Assorted Posts, 5 each (ø 1.1 mm, 1.3 mm, 1.5 mm)	5283