



hy-drop

happy
drop!

NEW
design

Smart and hygienic
dosage drop by drop

- New design
- One-hand operation
- No cross-contamination
- Made of high-tech polymer
- Excellent chemical resistance
- Autoclavable at 134 °C (273 °F)

POLYDENTIA SA
Swiss Products for Fine Dentistry



www.polydentia.com

With Hy-Drop save
liquids and time!



Pour the liquid into the
container (max 4ml).



Position the membrane.








Press gently on the
membrane. Hy-Drop releases
the liquid drop by drop.



Smart and hygienic liquid dispenser

ORDER FORM **Hy-Drop**

Please send me the following products:

Description	Quantity
 REF 5411/Y Hy-Drop yellow, 5 membranes	
 REF 5411/O Hy-Drop orange, 5 membranes	
 REF 5411/G Hy-Drop green, 5 membranes	
 REF 5411/V Hy-Drop violet, 5 membranes	
 REF 5411/B Hy-Drop blue, 5 membranes	
REF 5415 5 membranes, refill pack	

Name: _____

Address: _____

My dealer's name: _____

POLYDENTIA SA 
 Swiss Products for Fine Dentistry

DE Pat. 197 37 749
EP Pat. 1 006 852
US Pat. 6 267 270

6805 Mezzovico • Switzerland

Phone +41 (0)91 946 29 48

Fax +41 (0)91 946 32 03

www.polydentia.com

info@polydentia.com