

# Wood wedges

**Top Swiss quality, splinter-free interdental wedges**

- Anatomical shape, rounded tip
- Small size (XS, pinky) for paediatric dentistry
- Big size (XL) for periodontology
- Orange Wood Wedges - the most used size for your daily work

Square end  
for a secure grip

Rounded tip to  
prevent gingival trauma

Concave sides perfectly adapt  
along the entire cervical margin in the  
interdental space



Achievement of perfect contact points






Pinky 11 mm, the smallest wood wedge ever - ideal for paediatric dentistry













*Please send me the following products:*

Description and content	Type	REF	Q.ty
Wood Wedges, assorted, 100 pcs.	(XS, S, M, L)	5001	
Wood Wedges, assorted, 400 pcs.	(XS, S, M, L)	5040	
Wood Wedges, assorted, 600 pcs.	(XS, S, M, L)	5060	
Wood Wedges, assorted, 1000 pcs.	(XS, S, M, L)	5001+	

Complementary products	Type	REF	Q.ty
 Orange Wood Wedges, 200 pcs.	(SO)	5018	
 Wood Wedges, 100 pcs.	17 mm (XL)	5017	
 Wood Wedges, 1000 pcs.	17 mm (XL)	5017+	

Refills Wood Wedges	Type	REF	Q.ty
 Pinky Wood Wedges, 200 pcs.	11 mm (XS)	5011	
 Pinky Wood Wedges, 1000 pcs.	11 mm (XS)	5011+	
 Wood Wedges, 100 pcs.	12 mm (S)	5012	
 Wood Wedges, 1000 pcs.	12 mm (S)	5012+	
 Wood Wedges, 100 pcs.	13 mm (M)	5013	
 Wood Wedges, 1000 pcs.	13 mm (M)	5013+	
 Wood Wedges, 100 pcs.	15 mm (L)	5015	
 Wood Wedges, 1000 pcs.	15 mm (L)	5015+	

Address and signature:

My dealer is:

**Send your order by fax to:**  
**+41 91 946 32 03**

**or by e-mail to:**  
**sales@polydentia.com**

REV. 02 10.2014

